



INTRODUCTION TO BIOCONDUCTOR

Introducing ShortRead

Paula Andrea Martinez, PhD.
Data Scientist

Plant genomes

- *Arabidopsis thaliana* is a small flowering plant
- First plant to have its genome sequenced
- Genome size 135 megabase pairs (Mbp)





Sequencing companies





fastq vs fasta

fastq

```
@ unique sequence identifier  
raw sequence string  
  
+ optional id  
quality encoding per sequence letter
```

- fastq, fq

fasta

```
> unique sequence identifier  
raw sequence string
```

- fasta, fa, seq

fasta

```
library(ShortRead)

# read fasta
fasample <- readFasta(dirPath = "data/", pattern = "fasta")

# print fasample
print(fasample)
```

```
class: ShortRead
length: 500 reads; width: 50 cycles
```

```
# methods accessors
methods(class = "ShortRead")

# Write a ShortRead object
writeFasta(fasample, file = "data/sample.fasta")
```



fastq

```
library(ShortRead)

# read fastq
fqsample <- readFastq(dirPath = "data/", pattern = "fastq")

# print fqsample
fqsample
```

```
class: ShortReadQ
length: 500 reads; width: 50 cycles
```

```
# methods accessors
methods(class = "ShortReadQ")

# Write a ShortRead object
writeFastq(fqsample, file = "data/sample.fastq.gz")
```



fastq sample

```
library(ShortRead)

# set the seed to draw the same read sequences every time
set.seed(123)

# Subsample of 500 bases
sampler <- FastqSampler("data/SRR1971253.fastq", 500)

# save the yield of 500 read sequences
sample_small <- yield(sampler)

# Class ShortReadQ
class(sample_small)

# length 500 reads
length(sample_small)
```



INTRODUCTION TO BIOCONDUCTOR

You are ready!



INTRODUCTION TO BIOCONDUCTOR

Sequence quality

Paula Andrea Martinez, PhD.
Data Scientist



Quality scores - Phred table

| Quality value | Chance is wrong | Accuracy (%) |
|---------------|-----------------|--------------|
| 10 | 1 in 10 | 90 |
| 20 | 1 in 100 | 99 |
| 30 | 1 in 1 000 | 99.9 |
| 40 | 1 in 10 000 | 99.99 |
| 50 | 1 in 100 000 | 99.999 |



Encoding - Phred +33

```
# quality encoding
encoding(quality(fqsample))
```

Encoding characters and their scores

| | | | | | | | | | | | | | | | |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|---|----------|
| ! | " | # | \$ | % | & | ' | (|) | * | + | , | - | . | # | encoding |
| 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | # | score |
| / | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | : | ; | < | # | encoding |
| 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 | 27 | # | score |
| = | > | ? | @ | A | B | C | D | E | F | G | H | I | | # | encoding |
| 28 | 29 | 30 | 31 | 32 | 33 | 34 | 35 | 36 | 37 | 38 | 39 | 40 | | # | score |



fastq quality

```
library(ShortRead)
```

```
quality(fqsample)
```

```
class: FastqQuality  
A BStringSet instance
```

```
# Quality is represented with ASCII characters  
[1] 40 ?@@DDDDDDHDFDHE>AHFEGFIIEBGDBHH<3FEBEEEEEG  
[2] 40 BCCDFFFFHHHHHHJJJJJJJJJEHHGHIJJJJJJJJJJJ  
[3] 40 BCCFFFFFFHFHHHJJJJJJJIIJJIIIIIGIIJJIJGIJII  
[4] 40 CCCFFFFFFHHHHHHJJJJJJJJJJJJJJJJJJJJJJJJJJJJJJ
```



Exploring quality encoding and scores

```
library(ShortRead)
```

```
sread(fqsample)[1]
```

```
[1]      50 GTCCCATTTACCTCTGACTCTTTTGATGCTGCAATTGCTGCTCATATACT
```

```
# Quality is represented with ASCII characters  
quality(fqsample)[1]
```

```
[1]      50 ?@@DDDDDDHDFDHE>AHFEGFIIEBGDBHH<3FEBEEEEGGIGIIGHGHC
```

```
## PhredQuality instance  
pq <- PhredQuality(quality(fqsample))
```

```
# transform encoding into scores  
qs <- as(pq, "IntegerList")  
qs # print scores
```

```
30 31 31 35 35 35 35 35 39 35 37 35 39 36 29 32 39 37 36 38  
37 40 40 36 33 38 35 33 39 39 27 18 37 36 33 36 36 36 36 38  
38 40 38 40 40 38 39 38 39 34
```

Quality assessment

```
library(ShortRead)

# Quality assessment
qaSummary <- qa(fqsample, lane = 1)    # optional lane
```

```
# class: ShortReadQQA(10)
# Names accessible with the quality assessment summary
names(qaSummary)
```

```
[1] "readCounts"          "baseCalls"          "readQualityScore"
[4] "baseQuality"         "alignQuality"       "frequentSequences"
[7] "sequenceDistribution" "perCycle"           "perTile"
[10] "adapterContamination"
# QA elements are accessed with qa[["name"]]
```

```
# Get a HTML report
browseURL(report(qaSummary))
```



Alphabet by cycle

```
library(ShortRead)

# sequences alphabet
alphabet(sread(fullSample))
# [1] A,C,G,T,M,R,W,S,Y,K,V,H,D,B,N,-,+,.

abc <- alphabetByCycle(sread(fullSample))

# each observation is a letter and each variable is a cycle
# first, select the four first rows nucleotides A, C, G, T
# then, transpose
nucByCycle <- t(abc[1:4,])

nucByCycle <- nucByCycle %>%
  as.tibble() %>% # convert to tibble
  mutate(cycle = 1:50) # add cycle numbers

nucByCycle
```

| | A | C | G | T | cycle |
|-------|-------|-------|-------|---|-------|
| 16839 | 16335 | 16740 | 10878 | 1 | |
| 13056 | 13327 | 12064 | 22389 | 2 | |
| 13666 | 15617 | 13198 | 18355 | 3 | |
| 14723 | 15439 | 14239 | 16435 | 4 | |



INTRODUCTION TO BIOCONDUCTOR

Are you excited?



INTRODUCTION TO BIOCONDUCTOR

Match and filter

Paula Andrea Martinez, PhD.
Data Scientist



Duplicate sequences

- Biological sequence duplicates occur in nature
- Amplification from the steps in library preparation (PCR)
- Sequencing the sample more than once

Remove duplicates or at least mark them

- Whole genome sequencing or exome sequencing

Mark duplicates using a threshold

- RNA-seq and ChIP-seq



sruplicated

```
library(ShortRead)
```

```
# Counting duplicates TRUE is the number of duplicates  
table(sruplicated(dfqsample))
```

```
FALSE  TRUE  
500    500
```

```
# Cleaning reads from duplicates x[fun(x)]  
cleanReads <- mydReads[sruplicated(mydReads) == FALSE]
```

```
# Counting duplicates  
table(sruplicated(cleanReads))
```

```
FALSE  
500
```



Creating your own filters

`srFilter` to filter based on a condition `x[fun(x)]`

Filter example

```
library(SHORTREAD)

# Use a custom filter to remove reads from fqsample
# This filter to remove reads shorter than a min number of bases
readWidthCutoff <- srFilter(function(x) {width(x) >= minWidth},
                             name = "MinWidth")

minWidth <- 51

fqsample[readWidthCutoff(fqsample)]
```



nFilter

```
library(ShortRead)

# save your filter, .name is optional
myFilter <- nFilter(threshold = 10, .name = "cleanNFilter")

# use the filter at reading point
filtered <- readFastq(dirPath = "data",
                     pattern = ".fastq",
                     filter = myFilter)

# you will retrieve only those reads that have a maximum of 10 N's
filtered
```

idFilter and polynFilter

```
library(ShortRead)

#id filter example
myFilterID <- idFilter(regex = ":3:1")
# will return only those ids that contain the regular expression
# optional parameters are .name, fixed and exclude

# use the filter at reading point
filtered <- readFastq(dirPath = "data", pattern = ".fastq",
                     filter = myFilterID)
```

```
# filter to remove poly-A regions
myFilterPolyA <- polynFilter(threshold = 10, nuc = c("A"))
# will return the sequences that have a maximum number of 10 consecutive A's

# use the filter for subsetting
filtered[myFilterPolyA(filtered)]
```



INTRODUCTION TO BIOCONDUCTOR

Let's practice using filters!



INTRODUCTION TO BIOCONDUCTOR

Multiple and parallel sequence quality assessment

Paula Andrea Martinez, PhD.
Data Scientist

Rqc

```
library(Rqc)
```

- Uses Bioconductor packages that you have already used:
 - Biostrings, IRanges, methods, S4vectors
- New packages to discover in the following Bioconductor courses:
 - Rsamtools, GenomicAlignments, GenomicFiles, BiocParallel
- CRAN packages:
 - Knitr, dplyr, markdown, ggplot2, digest, shiny and Rccp



rqcQA

```
library(Rqc)

files <- # get the full path of the files you want to assess

qaRqc <- rqcQA(files)

# exploring qaRqc
class(qaRqc) # "list"
names(qaRqc) # name of the input files

# for each file
qaRqc[1]
# the class of the results is RqcResultSet
```

rqcQA arguments

```
library(Rqc)

# get the path of the files you want to assess
files <- "data/seq1.fq" "data/seq2.fq" "data/seq3.fq" "data/seq4.fq"

qaRqc <- rqcQA(files, workers = 4))

# sample of sequences
set.seed(1111)

qaRqc_sample <- rqcQA(files, workers = 4, sample = TRUE, n = 500))

# paired-end files
pfiles <- "data/seq_11.fq" "data/seq1_2.fq" "data/seq2_1.fq" "data/seq2_2.fq"

qaRqc_paired <- rqcQA(pfiles, workers = 4, pair = c(1, 1, 2, 2)))
```



rqcReport and rqcResultSet

```
# create a report
reportFile <- rqcReport(qaRqc, templateFile = "myReport.Rmd")

browseURL(reportFile)

#The class of qaRqc is rqcResultSet
methods(class = "RqcResultSet")
```


perFileInformation

```
qaRqc <- rqcQA(files, workers = 4))  
perFileInformation(qaRqc)
```

| filename | pair | format | group | reads | total.reads | path |
|------------------|------|--------|-------|-------|-------------|--------|
| SRR7760274.fastq | 1 | FASTQ | None | 1e+06 | 2404795 | ./data |
| SRR7760275.fastq | 2 | FASTQ | None | 1e+06 | 1508139 | ./data |
| SRR7760276.fastq | 3 | FASTQ | None | 1e+06 | 1950463 | ./data |
| SRR7760277.fastq | 4 | FASTQ | None | 1e+06 | 2629588 | ./data |

Plot functions

| rqc Plot functions | rqc Plot functions |
|--|--|
| <code>rqcCycleAverageQualityPcaPlot()</code> | <code>rqcGroupCycleAverageQualityPlot()</code> |
| <code>rqcCycleAverageQualityPlot()</code> | <code>rqcReadQualityBoxPlot()</code> |
| <code>rqcCycleBaseCallsLinePlot()</code> | <code>rqcReadQualityPlot()</code> |
| <code>rqcCycleBaseCallsPlot()</code> | <code>rqcReadWidthPlot()</code> |
| <code>rqcCycleGCPlot()</code> | <code>rqcReadFrequencyPlot()</code> |
| <code>rqcCycleQualityBoxPlot()</code> | <code>rqcCycleQualityPlot()</code> |





**KEEP
CALM**

AND use

**parallel quality
assessment**



INTRODUCTION TO BIOCONDUCTOR

You are ready!



INTRODUCTION TO BIOCONDUCTOR

Congratulations!

Paula Andrea Martinez, PhD.
Data Scientist

You learned...

- Install packages from Bioconductor by using the `BiocInstaller` package.
- Techniques for reading, manipulating and filtering raw genomic data using `BioStrings`, `GenomicRanges` **and** `ShortRead`.
- To work with `BSgenome` **and** `TxDb` built-in datasets. Then used these to identify patterns by using matching functions.
- Check the quality of sequence files using `ShortRead` **and** `Rqc`.

You explored





INTRODUCTION TO BIOCONDUCTOR

Keep learning!